


Custom Abutment Emergence Profile Options



NO Tissue Displacement

- **Open gingival embrasure**
- **Follow healing abutment emergence angle with no pressure on tissue**

SLIGHT Tissue Displacement

- **Natural gingival embrasure**
- **Follow healing abutment emergence angle with minimal pressure on tissue**

HEAVY Tissue Displacement

- **Closed gingival embrasure**
- **Follow healing abutment emergence angle with heavy pressure on tissue**

Margin Depth (Subgingival):

Anterior: Mesial 0.5mm Distal 0.5mm Lingual 0.3mm Buccal 1.0mm
Posterior: Mesial 0.5mm Distal 0.5mm Lingual 0.5mm Buccal 0.5mm

Default Custom Abutment Type:

Anterior: Titanium Gold Hue Titanium/Zirconia Hybrid
Posterior: Titanium Gold Hue Titanium/Zirconia Hybrid